



DENTAL PRODUCTS REPORT™

DISCOVERY. ADOPTION. SUCCESS.

dentalproductsreport.com

DENTSPLY Midwest



A CLASSIC RELAUNCHED

Highspeed handpiece without the stalling or constant adjustments.

MORE ON PAGE 3

*** COVER PRODUCTS**
are chosen based on innovation or their perceived impact on improving patient care and practice efficiency.

Cover Products and other
EDITOR'S CHOICE
products throughout this issue are chosen based on innovation or their perceived impact on improving patient care and practice efficiency.



HANDPIECE WITH SPEED AND CONTROL

The lightweight, ergonomic **Midwest® Stylus™ ATC** offers Speed-Sensing Intelligence (SSI) that constantly maintains bur speed according to load and automatically optimizes power to virtually eliminate stalling and the need to make frequent adjustments. Superior Turbine Suspension (STS) allows the handpiece to operate at speeds of 330,000 rpm without noticeable bur deflection and chattering for optimal control, precision and transference of power to the bur.

*Source: BPA Publisher Statement, December 2008.